

KD-Validated Anti-UHRF2 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG2184

Product Overview

Name	KD-Validated Anti-UHRF2 Mouse Monoclonal Antibody
Accession(Primary)	Q96PU4
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human UHRF2
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 90kDa, observed, 95 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:400-1:2,000Flow Cytometry (FCM): 1:2,000


Target information (Q96PU4)

Gene Name	UHRF2
About Gene Symbol	UHRF2 (Ubiquitin Like With PHD And Ring Finger Domains 2) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellula



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-UHRF2 antibody (Cat#ABG2184). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-UHRF2 antibody (Cat#ABG2184, 1:2,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-UHRF2 antibody (Cat#ABG2184). UHRF2 expression in wild-type (WT) and UHRF2 knockdown (KD) HSHC cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-UHRF2 antibody (Cat#ABG2184, 1:2,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.