

KD-Validated Anti-RNF31 Mouse Oligoclonal Antibody

Mouse oligoclonal antibody


ABG2188

Product Overview


Name	KD-Validated Anti-RNF31 Mouse Oligoclonal Antibody
Accession(Primary)	Q96EP0
Host Species	Mouse
Clonality	Oligoclonal
Isotype	Mouse IgG
Immunogen	Recombinant protein of human RNF31
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 120 kDa, observed, 120 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500Flow Cytometry (FCM): 1:2,000Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q96EP0)


Gene Name	RNF31
About Gene Symbol	RNF31 (Ring Finger Protein 31) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC



Immunocytochemical staining of HT-1080 cells with anti-RNF31 antibody (Cat#ABG2188, 1:1,000). Nuclei were stained blue with DAPI; RNF31 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Low. Scale bar, 20 μ m.