

KD-Validated Anti-Legumain Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2191

Product Overview


Name	KD-Validated Anti-Legumain Rabbit Monoclonal Antibody
Accession(Primary)	Q99538
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Legumain
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 49 kDa, observed, 36,57 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q99538)


Gene Name	LGMN
About Gene Symbol	LGMN (Legumain) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2025 Genuin Biotechnologies LLC



ht ©2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC