

KD-Validated Anti-REEP2 Mouse Monoclonal Antibody

Mouse monoclonal Antibody


ABG2193

Product Overview


Name	KD-Validated Anti-REEP2 Mouse Monoclonal Antibody
Accession(Primary)	Q9BRK0
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Immunogen	Recombinant protein of human REEP2
Species Reactivity	Human
Calculated MW	Predicted, 28 kDa, observed, 28 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (Q9BRK0)

Gene Name	REEP2
About Gene Symbol	REEP2 (Receptor Accessory Protein 2) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2025 Genuin Biotechnologies LLC



Copyright © 2025 Genuin Biotechnologies LLC