

KD-Validated Anti-CENPA Rabbit Monoclonal Antibody

Rabbit monoclonal Antibody


ABG2196

Product Overview


Name	KD-Validated Anti-CENPA Rabbit Monoclonal Antibody
Accession(Primary)	P49450
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CENPA
Species Reactivity	Human
Calculated MW	Predicted, 16 kDa, observed, 18 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P49450)

Gene Name	CENPA
About Gene Symbol	CENPA (Centromere Protein A) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics and adaptive immunity. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC