

KD-Validated Anti-COX5B Rabbit Monoclonal Antibody

Rabbit monoclonal Antibody


ABG2200

Product Overview

Name	KD-Validated Anti-COX5B Rabbit Monoclonal Antibody
Accession(Primary)	P10606
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human COX5B
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 14 kDa, observed, 14 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000Flow Cytometry (FCM): 1:1,000


Target information (P10606)

Gene Name	COX5B
About Gene Symbol	COX5B (Cytochrome C Oxidase Subunit 5B) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identi




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-COX5B antibody (Cat#ABG2200). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-COX5B antibody (Cat#ABG2200, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-COX5B antibody (Cat#ABG2200). COX5B expression in wild-type (WT) and COX5B shRNA knockdown (KD) HeLa cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-COX5B antibody (Cat#ABG2200, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of COX5B expression in H9c2 cells using anti-COX5B antibody (Cat#ABG2200, 1:2,000). Green, isotype control; red, COX5B.