

KD-Validated Anti-EFEMP1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2210

Product Overview

Name	KD-Validated Anti-EFEMP1 Rabbit Monoclonal Antibody
Accession(Primary)	Q12805
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	Recombinant protein of human EFEMP1/Fibulin-3
Species Reactivity	Human
Calculated MW	Predicted, 55 kDa; observed, 55-65 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000


Target information (Q12805)

Gene Name	EFEMP1
About Gene Symbol	EFEMP1 (EGF Containing Fibulin Extracellular Matrix Protein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-EFEMP1 antibody (Cat#ABG2210). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-EFEMP1 antibody (Cat#ABG2210, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-EFEMP1 antibody (Cat#ABG2210). EFEMP1 expression in wild-type (WT) and EFEMP1 shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-EFEMP1 antibody (Cat#ABG2210, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of EFEMP1 expression in HeLa cells using anti-EFEMP1 antibody (Cat#ABG2210, 1:2,000). Green, isotype control; red, EFEMP1.