

KD-Validated Anti-ZC3HC1 Mouse Oligoclonal Antibody

Mouse oligoclonal antibody


ABG2228

Product Overview


Name	KD-Validated Anti-ZC3HC1 Mouse Oligoclonal Antibody
Accession(Primary)	Q86WB0
Host Species	Mouse
Clonality	Oligoclonal
Isotype	Mouse IgG
Immunogen	Recombinant protein of human ZC3HC1
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 55 kDa; observed, 55 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500 Flow Cytometry (FCM): 1:2,000

Target information (Q86WB0)


Gene Name	ZC3HC1
About Gene Symbol	ZC3HC1 (Zinc Finger C3HC-Type Containing 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC