

KD-Validated Anti-TFAP2C Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2229

Product Overview

Name	KD-Validated Anti-TFAP2C Rabbit Monoclonal Antibody
Accession(Primary)	Q92754
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	Recombinant protein of human TFAP2C
Species Reactivity	Human
Calculated MW	Predicted, 49 kDa; observed, 49 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (Q92754)

Gene Name	TFAP2C
About Gene Symbol	TFAP2C (Transcription Factor AP-2 Gamma) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-TFAP2C antibody (Cat#ABG2229). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-TFAP2C antibody (Cat#ABG2229, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-TFAP2C antibody (Cat#ABG2229). TFAP2C expression in wild-type (WT) and TFAP2C shRNA knockdown (KD) HeLa cells with 30 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-TFAP2C antibody (Cat#ABG2229, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.