

KD-Validated Anti-Malic Enzyme 2 Rabbit Monoclonal

Rabbit monoclonal antibody


ABG2231

Product Overview

Name	KD-Validated Anti-Malic Enzyme 2 Rabbit Monoclonal Antibody
Accession(Primary)	P23368
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human ME2
Calculated MW	Predicted, 65 kDa; observed, 58 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (P23368)

Gene Name	ME2
About Gene Symbol	



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-malic enzyme 2 antibody (Cat#ABG2231). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-malic enzyme 2 antibody (Cat#ABG2231, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-malic enzyme 2 antibody (Cat#ABG2231). Malic enzyme 2 expression in wild-type (WT) and malic enzyme 2 (ME2) shRNA knockdown (KD) HeLa cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-malic enzyme 2 antibody (Cat#ABG2231, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.