

KD-Validated Anti-UBE2F Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2233

Product Overview

Name	KD-Validated Anti-UBE2F Rabbit Monoclonal Antibody
Accession(Primary)	Q969M7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	Recombinant protein of human UBE2F
Species Reactivity	Human
Calculated MW	Predicted, 21 kDa; observed, 20 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (Q969M7)

Gene Name	UBE2F
About Gene Symbol	UBE2F (Ubiquitin Conjugating Enzyme E2 F (Putative)) is a metabolic enzyme involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-UBE2F antibody (Cat#ABG2233). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-UBE2F antibody (Cat#ABG2233, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-UBE2F antibody (Cat#ABG2233). UBE2F expression in wild-type (WT) and UBE2F shRNA knockdown (KD) HT-1080 cells with 30 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-UBE2F antibody (Cat#ABG2233, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.