

KD-Validated Anti-Cullin 5 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2249

Product Overview

Name	KD-Validated Anti-Cullin 5 Rabbit Monoclonal Antibody
Accession(Primary)	Q93034
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	Recombinant protein of human Cullin 5
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 91 kDa; observed, 75 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (Q93034)

Gene Name	CUL5
About Gene Symbol	CUL5 (Cullin 5) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-cullin 5 antibody (Cat#ABG2249). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-cullin 5 antibody (Cat#ABG2249, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-cullin 5 antibody (Cat#ABG2249). Cullin 5 expression in wild-type (WT) and Cullin 5 (CUL5) shRNA knockdown (KD) 293T cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-cullin 5 antibody (Cat#ABG2249, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.