

KD-Validated Anti-Lactate Dehydrogenase C Rabbit

Monoclonal Antibody

Rabbit monoclonal antibody


ABG2257

Product Overview


Name	KD-Validated Anti-Lactate Dehydrogenase C Rabbit Monoclonal Antibody
Accession(Primary)	P07864
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Lactate Dehydrogenase C
Species Reactivity	Human
Calculated MW	Predicted, 36 kDa; observed, 31 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P07864)

Gene Name	LDHC
About Gene Symbol	LDHC (Lactate Dehydrogenase C) is a metabolic dehydrogenase involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC