

KD-Validated Anti-RAB10 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2263

Product Overview

Name	KD-Validated Anti-RAB10 Rabbit Monoclonal Antibody
Accession(Primary)	P61026
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	Recombinant protein of human RAB10
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 23 kDa; observed, 20 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (P61026)

Gene Name	RAB10
About Gene Symbol	RAB10 (RAB10, Member RAS Oncogene Family) is a small GTPase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-RAB10 antibody (Cat#ABG2263). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-RAB10 antibody (Cat#ABG2263, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-RAB10 antibody (Cat#ABG2263). RAB10 expression in wild-type (WT) and RAB10 shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-RAB10 antibody (Cat#ABG2263, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.