

KD-Validated Anti-TACSTD2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2271

Product Overview

Name	KD-Validated Anti-TACSTD2 Rabbit Monoclonal Antibody
Accession(Primary)	P09758
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TACSTD2
Species Reactivity	Human
Calculated MW	Predicted, 36 kDa; observed, 36-55 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (P09758)

Gene Name	TACSTD2
About Gene Symbol	TACSTD2 (Tumor Associated Calcium Signal Transducer 2) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify the



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-TACSTD2 antibody (Cat#ABG2271). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-TACSTD2 antibody (Cat#ABG2271, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-TACSTD2 antibody (Cat#ABG2271). TACSTD2 expression in wild-type (WT) and TACSTD2 shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-TACSTD2 antibody (Cat#ABG2271, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.