

KD-Validated Anti-TEAD4 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2273

Product Overview

Name	KD-Validated Anti-TEAD4 Rabbit Monoclonal Antibody
Accession(Primary)	Q15561
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	Recombinant protein of human TEAD4
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 48 kDa; observed, 54 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (Q15561)

Gene Name	TEAD4
About Gene Symbol	TEAD4 (TEA Domain Transcription Factor 4) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-TEAD4 antibody (Cat#ABG2273). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-TEAD4 antibody (Cat#ABG2273, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-TEAD4 antibody (Cat#ABG2273). TEAD4 expression in wild type (WT) and TEAD4 knockdown (KD) HepG2 cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-TEAD4 antibody (Cat#ABG2273, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.