

KD-Validated Anti-TXNDC9 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2276

Product Overview

Name	KD-Validated Anti-TXNDC9 Rabbit Monoclonal Antibody
Accession(Primary)	O14530
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TXNDC9
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 27 kDa; observed, 27 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (O14530)

Gene Name	TXNDC9
About Gene Symbol	TXNDC9 (Thioredoxin Domain Containing 9) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-TXNDC9 antibody (Cat#ABG2276). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-TXNDC9 antibody (Cat#ABG2276, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-TXNDC9 antibody (Cat#ABG2276). TXNDC9 expression in wild-type (WT) and TXNDC9 shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-TXNDC9 antibody (Cat#ABG2276, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.