

KD-Validated Anti-SNRPC Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2277

Product Overview

Name	KD-Validated Anti-SNRPC Rabbit Monoclonal Antibody
Accession(Primary)	P09234
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human SNRPC
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 17 kDa; observed, 19 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (P09234)

Gene Name	SNRPC
About Gene Symbol	SNRPC (Small Nuclear Ribonucleoprotein Polypeptide C) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-SNRPC antibody (Cat#ABG2277). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-SNRPC antibody (Cat#ABG2277, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-SNRPC antibody (Cat#ABG2277). SNRPC expression in wild-type (WT) and SNRPC shRNA knockdown (KD) HeLa cells with 30 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-SNRPC antibody (Cat#ABG2277, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.