

KD-Validated Anti-LLGL1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2281

Product Overview

Name	KD-Validated Anti-LLGL1 Rabbit Monoclonal Antibody
Accession(Primary)	Q15334
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	Recombinant protein of human LLGL1
Species Reactivity	Human
Calculated MW	Predicted, 115 kDa; observed, 130 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (Q15334)

Gene Name	LLGL1
About Gene Symbol	LLGL1 (LLGL Scribble Cell Polarity Complex Component 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-LLGL1 antibody (Cat#ABG2281). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-LLGL1 antibody (Cat#ABG2281, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-LLGL1 antibody (Cat#ABG2281). LLGL1 expression in wild-type (WT) and LLGL1 shRNA knockdown (KD) HeLa cells with 30 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-LLGL1 antibody (Cat#ABG2281, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.