

KD-Validated Anti-TPX2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2286

Product Overview

Name	KD-Validated Anti-TPX2 Rabbit Monoclonal Antibody
Accession(Primary)	Q9ULW0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TPX2
Species Reactivity	Human
Calculated MW	Predicted, 86 kDa; observed, 98 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (Q9ULW0)

Gene Name	TPX2
About Gene Symbol	TPX2 (TPX2 Microtubule Nucleation Factor) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targ



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-TPX2 antibody (Cat#ABG2286). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-TPX2 antibody (Cat#ABG2286, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-TPX2 antibody (Cat#ABG2286). TPX2 expression in wild-type (WT) and TPX2 shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-TPX2 antibody (Cat#ABG2286, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.