

KD-Validated Anti-active + pro caspase 3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG2292

Product Overview

Name	KD-Validated Anti-active + pro caspase 3 Rabbit Monoclonal Antibody
Accession(Primary)	P42574
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human active + pro Caspase 3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 32 kDa , observed, 32 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P42574)

Gene Name	CASP3
About Gene Symbol	CASP3 (Caspase 3) is a cysteine protease (caspase) involved in apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

