

KD-Validated Anti-Ephrin Receptor B4 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG2299

Product Overview

Name	KD-Validated Anti-Ephrin Receptor B4 Rabbit Monoclonal Antibody
Accession(Primary)	P54760
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Eph receptor B4/HTK
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted , 108 kDa , observed , 135 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (P54760)

Gene Name	EPHB4
About Gene Symbol	EPHB4 (EPH Receptor B4) is a receptor tyrosine kinase involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

