

KD-Validated Anti-Albumin Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2324

Product Overview


Name	KD-Validated Anti-Albumin Rabbit Monoclonal Antibody
Accession(Primary)	P02768
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from Albumin
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 69 kDa ; observed, 69 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P02768)

Gene Name	ALB
About Gene Symbol	ALB (Albumin) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-Albumin antibody (Cat#ABG2324). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Albumin antibody (Cat#ABG2324, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-Albumin antibody (Cat#ABG2324). Albumin expression in wild type (WT) and Albumin shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-Albumin antibody (Cat#ABG2324, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.