

KD-Validated Anti-PRAS40 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG2325

Product Overview

Name	KD-Validated Anti-PRAS40 Rabbit Monoclonal Antibody
Accession(Primary)	Q96B36
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PRAS40
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 27 kDa; Observed, 40 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:2,500 Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q96B36)

Gene Name	AKT1S1
About Gene Symbol	AKT1S1 (AKT1 Substrate 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

