

KD-Validated Anti-CD163 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2326

Product Overview


Name	KD-Validated Anti-CD163 Rabbit Monoclonal Antibody
Accession(Primary)	Q86VB7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CD163
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 125 kDa; Observed, 125 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (Q86VB7)


Gene Name	CD163
About Gene Symbol	CD163 (CD163 Molecule) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-CD163 antibody (Cat#ABG2326). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-CD163 antibody (Cat#ABG2326, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-CD163 antibody (Cat#ABG2326). CD163 expression in wild type (WT) and CD163 shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-CD163 antibody (Cat#ABG2326, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Immunocytochemical staining of HeLa cells with CD163 antibody (Cat#ABG2326, 1:1,000). Nuclei were stained blue with DAPI; CD163 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 µm.