

KD-Validated Anti-Ferritin heavy chain 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2339

Product Overview


Name	KD-Validated Anti-Ferritin heavy chain 1 Rabbit Monoclonal Antibody
Accession(Primary)	P02794
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Ferritin
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 21 kDa; Observed, 15 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Flow cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P02794)


Gene Name	FTH1
About Gene Symbol	FTH1 (Ferritin Heavy Chain 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with Ferritin heavy chain 1 antibody (Cat#ABG2339, 1:1,000). Nuclei were stained blue with DAPI; Ferritin heavy chain 1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.