

KD-Validated Anti-Phospho-GSK3 beta (Ser9) Rabbit Monoclonal Antibody


ABG2341

Product Overview


Name	KD-Validated Anti-Phospho-GSK3 beta (Ser9) Rabbit Monoclonal Antibody
Accession(Primary)	P49841
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-GSK3 beta (Ser9)
Species Reactivity	Human
Calculated MW	Predicted, 4 7 kDa ; Observed, 4 6 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:10,000 Immunocytochemistry (IC): 1:1,000

Target information (P49841)


Gene Name	GSK3B
About Gene Symbol	GSK3B (Glycogen Synthase Kinase 3 Beta) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC