

KD-Validated Anti-Angio associated migratory cell protein

Rabbit monoclonal antibody


ABG2346

Product Overview

Name	KD-Validated Anti-Angio associated migratory cell protein Rabbit Monoclonal Antibody
Accession(Primary)	Q13685
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human AAMP
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 47 kDa ; Observed, 52 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000


Target information (Q13685)

Gene Name	AAMP
About Gene Symbol	AAMP (Angio Associated Migratory Cell Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Angio associated migratory cell protein antibody (Cat#ABG2346). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Angio associated migratory cell protein antibody (Cat#ABG2346, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Angio associated migratory cell protein antibody (Cat#ABG2346). Angio associated migratory cell protein expression in wild type (WT) and angio associated migratory cell protein shRNA knockdown (KD) HeLa cells with 30 μ g of total cell lysates. β -Tubulin serves as a loading control. The blot was incubated with anti-Angio associated migratory cell protein antibody (Cat#ABG2346, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.