

KD-Validated Anti-Phospho-PKR (T451) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2357

Product Overview


Name	KD-Validated Anti-Phospho-PKR (T451) Rabbit Monoclonal Antibody
Accession(Primary)	P19525
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-PKR (T451)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 62 kDa ; Observed, 68 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (P19525)

Gene Name	EIF2AK2
About Gene Symbol	EIF2AK2 (eIF-2A protein kinase 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC

