

KD-Validated Anti-EHMT2/G9A Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2359

Product Overview


Name	KD-Validated Anti-EHMT2/G9A Rabbit Monoclonal Antibody
Accession(Primary)	Q96KQ7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human EHMT2/G9A
Species Reactivity	Human
Calculated MW	Predicted, 132 kDa; Observed, 160-170 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:2,000-:10,000 Flow cytometry (FCM): 1:2,000

Target information (Q96KQ7)


Gene Name	EHMT2
About Gene Symbol	EHMT2 (Euchromatic Histone Lysine Methyltransferase 2) is a histone methyltransferase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC