

KD-Validated Anti-ERCC1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG2374

Product Overview

Name	KD-Validated Anti-ERCC1 Rabbit Monoclonal Antibody
Accession(Primary)	P07992
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human ERCC1
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 33 kDa ; Observed, 36 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (P07992)

Gene Name	ERCC1
About Gene Symbol	ERCC1 (ERCC Excision Repair 1, Endonuclease Non-Catalytic Subunit) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

