

KD-Validated Anti-CRKII Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2377

Product Overview


Name	KD-Validated Anti-CRKII Rabbit Monoclonal Antibody
Accession(Primary)	P46108
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CRKII
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 34 kDa; Observed, 35 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Flow cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P46108)


Gene Name	CRK
About Gene Symbol	CRK (CRK Proto-Oncogene, Adaptor Protein) is a signaling adaptor protein involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with CRKII antibody (Cat#ABG2377, 1:1,000). Nuclei were stained blue with DAPI; CRKII was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.