

KD-Validated Anti-Stat3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody





ABG2389


Product Overview

Name	KD-Validated Anti-Stat3 Rabbit Monoclonal Antibody
Accession(Primary)	P40763
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human beta STAT3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 88 kDa ; Observed, 88 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Flow cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P40763)

Gene Name	STAT3
About Gene Symbol	STAT3 (Signal Transducer And Activator Of Transcription 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with anti-Stat3Antibody antibody (Cat#ABG2389, 1:1,000). Nuclei were stained blue with DAPI; Stat3Antibody was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.