

KD-Validated Anti-C-terminal Src kinase Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2393

Product Overview


Name	KD-Validated Anti-C-terminal Src kinase Rabbit Monoclonal Antibody
Accession(Primary)	P41240
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CSK
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 51 kDa; Observed, 48 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (P41240)


Gene Name	CSK
About Gene Symbol	CSK (C-Terminal Src Kinase) is a receptor tyrosine kinase involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genin Biotechnologies LLC



Copyright ©2024 Genin Biotechnologies LLC



Immunocytochemical staining of HT-1080 cells with C-terminal Src kinase antibody (Cat#ABG2393, 1:1,000). Nuclei were stained blue with DAPI; C-terminal Src kinase was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 µm.