

KO Validated Anti-BCL2 Like 11 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2411

Product Overview

Name	KO Validated Anti-BCL2 Like 11 Rabbit Monoclonal Antibody
Accession(Primary)	O43521
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Bim
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 22 kDa, observed, 22 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000


Target information (O43521)

Gene Name	BCL2L11
About Gene Symbol	BCL2L11 (BCL2 Like 11) is a structural cytoskeletal protein involved in apoptosis. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-BCL2 like 11 antibody (Cat#ABG2411). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-BCL2 like 11 antibody (Cat#ABG2411, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-BCL2 like 11 antibody (Cat#ABG2411). BCL2 like 11 expression in wild type (WT) and BCL2 like 11 (BCL2L11) knockout (KO) HSHC cells with 20 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-BCL2 like 11 antibody (Cat#ABG2411, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of BCL2 like 11 expression in HepG2 cells using anti-BCL2 like 11 antibody (Cat# ABG2411, 1:2,000). Green, isotype control; red, BCL2 like 11.