

KO-Validated Anti-HADH Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG2425

Product Overview

Name	KO-Validated Anti-HADH Mouse Monoclonal Antibody
Accession(Primary)	Q16836
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Immunogen	Recombinant protein of human HADH
Species Reactivity	Human
Calculated MW	Predicted, 34 kDa; observed, 30 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:400-1:2,000


Target information (Q16836)

Gene Name	HADH
About Gene Symbol	HADH (Hydroxyacyl-CoA Dehydrogenase) is a metabolic dehydrogenase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-HADH antibody (Cat#ABG2425). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-HADH antibody (Cat#ABG2425, 1:2,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-HADH antibody (Cat#ABG2425). HADH expression in wild type (WT) and HADH knockout (KO) HeLa cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-HADH antibody (Cat#ABG2425, 1:2,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.