

KO-Validated Anti-ETS Transcription Factor ELK1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2428

Product Overview

Name	KO-Validated Anti-ETS Transcription Factor ELK1 Rabbit Monoclonal Antibody
Accession(Primary)	P19419
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human ELK1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 45 kDa; observed, 62 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000Flow Cytometry (FCM): 1:2,000Immunocytochemistry (IC): 1:100-1:1,000


Target information (P19419)

Gene Name	ELK1
About Gene Symbol	ELK1 (ETS Transcription Factor ELK1) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-ELK1 antibody (Cat#ABG2428). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-ELK1 antibody (Cat#ABG2428, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-ELK1 antibody (Cat#ABG2428). ELK1 expression in wild type (WT) and ELK1 knockout (KO) HSHC cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-ELK1 antibody (Cat#ABG2428, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of ELK1 expression in HepG2 cells using anti-ELK1 antibody (Cat#ABG2428, 1:2,000). Green, isotype control; red, ELK1.



Immunocytochemical staining of HepG2 cells with anti-ELK1 antibody (Cat#ABG2428, 1:1,000).

Nuclei were stained blue with DAPI; ELK1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.