

KO-Validated Anti-Ribosomal Protein SA Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2429

Product Overview


Name	KO-Validated Anti-Ribosomal Protein SA Rabbit Monoclonal Antibody
Accession(Primary)	P08865
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human 67kDa Laminin Receptor
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 33 kDa; observed, 42 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P08865)


Gene Name	RPSA
About Gene Symbol	RPSA (Ribosomal Protein SA) is an extracellular matrix laminin protein involved in translational regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC




Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC

Ribosomal protein SA-Alexa Fluor® 647



Immunocytochemical staining of HepG2 cells with anti-Ribosomal protein SA antibody (Cat#ABG2429, 1:1,000). Nuclei were stained blue with DAPI; Ribosomal protein SA was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.