

KO-Validated Anti-Phospho-Smad3 (S423 + S425) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2433

Product Overview


Name	KO-Validated Anti-Phospho-Smad3 (S423 + S425) Rabbit Monoclonal Antibody
Accession(Primary)	P84022
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-Smad3 (S423 + S425)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 61 kDa; observed, 60 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P84022)


Gene Name	SMAD3
About Gene Symbol	SMAD3 (SMAD Family Member 3) is a transcription factor involved in TGF- β signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC