

KO-Validated Anti-Phospho-YAP1 (S127) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2436

Product Overview


Name	KO-Validated Anti-Phospho-YAP1 (S127) Rabbit Monoclonal Antibody
Accession(Primary)	P46937
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human YAP1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 54 kDa; observed, 65 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:4,000-1:20,000Flow Cytometry (FCM): 1:2,000Immunocytochemistry (IC): 1:100-1:1,000

Target information (P46937)


Gene Name	YAP1
About Gene Symbol	YAP1 (Yes1 Associated Transcriptional Regulator) is a transcriptional coactivator involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC




Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC

Phospho-YAP1 (S127)-Alexa Fluor® 647



Immunocytochemical staining of HepG2 cells with anti-Phospho-YAP1 (S127) antibody (Cat#ABG2436, 1:1,000). Nuclei were stained blue with DAPI; Phospho-YAP1 (S127) was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.