

KO-Validated Anti-FTO Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG2440

Product Overview


Name	KO-Validated Anti-FTO Mouse Monoclonal Antibody
Accession(Primary)	Q9C0B1
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human FTO
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 58 kDa; observed, 58 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:500-1:2,500 Flow Cytometry (FCM): 1:2,000

Target information (Q9C0B1)


Gene Name	FTO
About Gene Symbol	FTO (FTO Alpha-Ketoglutarate Dependent Dioxygenase) is a metabolic enzyme involved in translational regulation. It is an important research target in metabolic disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC