

## EGFP Antibody

*Enhanced Green Fluorescent Protein, Mouse Anti*

ABM0218

### Product Overview

|                    |                                                                                                                                                                                                  |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>        | EGFP Antibody                                                                                                                                                                                    |
| <b>Immunogen</b>   | Anti-human EGFR mAb, is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice immunized with recombinant human EGFR amino acids 424-605 purified from E. coli. |
| <b>Description</b> | Enhanced Green Fluorescent Protein, Mouse Anti                                                                                                                                                   |