

# HA-tag Antibody

*HA-tag peptide, Mouse Antibody*

ABM0304

## Product Overview

|                    |                                                                                                                                                               |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>        | HA-tag Antibody                                                                                                                                               |
| <b>Immunogen</b>   | Anti-human HA-tag Peptide mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a synthetic peptide. |
| <b>Description</b> | HA-tag peptide, Mouse Antibody                                                                                                                                |