

HSP90B1 P94B9AT Antibody

Heat Shock Protein 90kDa Beta (GRP94) Member 1 , Mouse Anti Human

ABM0345

Product Overview

Name	HSP90B1 P94B9AT Antibody
Immunogen	Anti-human HSPA13 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human HSPA13 protein 23-471 amino acids purified from E. coli.
Description	Heat Shock Protein 90kDa Beta (GRP94) Member 1 , Mouse Anti Human