

SARS CoV-2 IgG S1

Recombinant Anti Human SARS CoV-2 IgG Spike S1, Monoclonal

ABM0606

Product Overview

Name	SARS CoV-2 IgG S1
Immunogen	The native monoclonal antibody was generated by sequencing peripheral blood lymphocytes of a patient exposed to the SARS-CoV. Mouse Anti Human SARS CoV-2 IgM Kappa Spike S1 Reformatted using the variable domain sequences of the original Human IgG1 format.
Description	Recombinant Anti Human SARS CoV-2 IgG Spike S1, Monoclonal