

BMP-7 Antibody

Bone Morphogenetic Protein-7 Polyclonal Rabbit Anti Human Antibody

ABP0004

Product Overview

Name	BMP-7 Antibody
Immunogen	Anti-human BMP8B mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human BMP8B protein 264-402 amino acids purified from E. coli.
Description	Bone Morphogenetic Protein-7 Polyclonal Rabbit Anti Human Antibody