

[KO] DAZAP1 Rabbit mAb

ABZ1022

Product Overview

Name	[KO] DAZAP1 Rabbit mAb
Accession(Primary)	Q96EP5
Host Species	Rabbit
Isotype	IgG
Immunogen	Recombinant protein of human DAZAP1
Species Reactivity	Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB: 1/500-1/1000; IP: 1/20

Target information (Q96EP5)

Gene Name	DAZAP1
About Gene Symbol	DAZAP1 (DAZ Associated Protein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.