

[KO] GRK6 Rabbit mAb

ABZ1050

Product Overview

Name	[KO] GRK6 Rabbit mAb
Accession(Primary)	P43250
Host Species	Rabbit
Isotype	IgG
Immunogen	A synthetic peptide of human GRK6
Calculated MW	Calculated MW: 66 kDa; Observed MW: 66 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (P43250)

Gene Name	GRK6
About Gene Symbol	GRK6 (G Protein-Coupled Receptor Kinase 6) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.