

[KO] N Cadherin Mouse mAb

ABZ1063

Product Overview

Name	[KO] N Cadherin Mouse mAb
Accession(Primary)	P19022
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human CDH2 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 100 kDa; Observed MW: 120 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200

Target information (P19022)

Gene Name	CDH2
About Gene Symbol	