

[KO] Cyclophilin F Mouse mAb

ABZ1065

Product Overview

Name	[KO] Cyclophilin F Mouse mAb
Accession(Primary)	P30405
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human PPIF produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P30405)

Gene Name	PPIF
About Gene Symbol	PPIF (Peptidylprolyl Isomerase F) is a peptidylprolyl isomerase (cyclophilin) involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.